

---

## PRODUCT SHEET



### Arctiko MUF40 Ultralow-Freezer (-80°)

The MUF 40 -80°C ultralow-freezer can be used both as a bench top and under counter model. It features a compact design and low noise level, making it ideal for operating inside the laboratory.

#### TECHNICAL DATA:

---

Volume: 35 l

Temperature: -60°C to -80°C

Dimensions exterior (WxDxH): 495 x 660 x 680 mm

Dimensions interior (WxDxH): 336 x 356 x 295 mm

Weight: 60 kg

#### SERVICE:



Calibration

#### SERVICE:



Temperature Mapping

#### SERVICE:



IQ/OQ

#### ANCILLARIES:



Temperatur Datalogger  
EBI 310-TPX220