
PRODUCT SHEET



Thermo HFU600TV large capacity Ultralow-Freezer (-86°)

The HFU600TV ultralow-freezer features a large 815 liter volume with a storage capacity for 600 standard cryoboxes.

TECHNICAL DATA:

Volume: 815 l

Temperature range: -50° C to -86° C

Exterior dimensions (WxDxH): 111 x 96 x 198 cm

Interior dimensions (WxDxH): 87 x 72 x 130 cm

Weight: 388 kg

SERVICE:



Calibration

SERVICE:



Temperature Mapping

SERVICE:



IQ/OQ

ANCILLARIES:



Temperatur Datalogger
EBI 310-TPX220