

**LIEBHERR**

# MKUv 1610 mit Schubladen-Ausstattung

Gross capacity, total	<b>142 l</b>
Product dimensions (H/W/D)	<b>82 / 60.1 / 61.8 cm</b>
Energy consumption in 24h	<b>0.747 kWh / 24h</b>
SmartMonitoring	<b>prepared</b>

> [Get help for service](#)



## General product information

Product group output channels

**Pharmacy refrigerator with Comfort controller**

GTIN

**9005382238158**

## Output and consumption

Energy consumption in 24h

**0.747 kWh / 24h**

Noise level

**47 dB**

Net capacity, total	<b>76 l</b>
Gross capacity, total	<b>142 l</b>
Refrigerant	<b>R 600a</b>
Voltage	<b>220-240 V ~</b>
Frequency	<b>50 Hz</b>
Connection rating	<b>1.0 A</b>

## Dimensions and weight

Product dimensions (H/W/D)	<b>82 / 60.1 / 61.8 cm</b>
Insulation	<b>40 - 40 mm</b>
Weight (with packaging)	<b>44.00 kg</b>
Weight (without packaging)	<b>42.00 kg</b>

## Control and functions

Control	<b>Electronic control system</b>
Control unit	<b>Buttons</b>
Malfunction: Warning signal	<b>optical and acoustical</b>
potential-free contact	<b>✓</b>
Power failure alarm	<b>immediately upon power failure for min. 12h</b>
Maximum fluctuation	<b>1.3 °C</b>
Gradient	<b>2.8 °C</b>
SmartMonitoring	<b>prepared</b>

## Design and features

Side wall material	<b>steel</b>
Colour	<b>white</b>
Door/Cover material	<b>steel</b>
Lock	<b>fitted</b>
Self-closing door	<b>✓</b>

Door hinges	<b>right interchangeable</b>
Material of interior containers	<b>Plastic white</b>
Type of castors	—
Door design	<b>HardLine</b>
Handle	<b>Ergonomic slimline handle with opening mechanism</b>

Please note: despite careful data management, we reserve the right to make technical changes; and discrepancies may occur between the image or text content and the original appliance.

---

**LIEBHERR**