
PRODUCT SHEET



Cardio M 12-channel ECG

The Cardio M is a reliable, easy-to-use 12-channel ECG which is one of our standard devices in clinical studies.

TECHNICAL DATA:

12-Channel Resting-ECG

Parameters: Heart rate, PR, QRS, QT/QTc, P-R-T axes

HRV-analysis: no

Clear print-out

Manual ECG recording or 10-sec ECG for interpretation

ECG diagram can easily be printed out several times for archiving purposes

Dimensions (WxDxH): 296 x 305 x 97 mm

Weight: 3,2 kg

SERVICE:



Annual Safety
Check

SERVICE:



Calibration

SERVICE:



reCAL Exchange
Program

SERVICE:



User Training

ANCILLARIES:



Adhesive Electrodes