
PRODUCT SHEET



Vitalograph Micro Spirometer

The Vitalograph Micro is a very compact spirometer with a simple, self explaining touchscreen operation. It fulfills the ATS/ERS Spirometer standard when using the 3 liter calibration pump.

Patient safety is ensured thanks to the single use Vitalograph BVF Bacterial-Viral-Filters.

TECHNICAL DATA:

Recorded values: VC, FVC, FEV1, FEV1/VK, PEF, FEF25-75, FEF25, FEF50, FEF75 and more than 30 other paramters (up to 8 can be chosen in parallel)

Calibration pump 3 liters ensures ATS/ERS spirometry standard

Dimensions (WxDxH): 83 x 93 x 32 mm

Weight: 250 grams

SERVICE:



Calibration

SERVICE:



reCAL Exchange
Program

ANCILLARIES:



Vitalograph BVF Bacterial-Viral-Filters
Vitalograph Disposable nose clips